

Provider 790 Fall Prevention Monitoring and Alerting

An integral part of your Patient Safety Fall Prevention Program

“81% of nurses spend at least 1-2 hours per shift checking beds to proactively ensure patient safety.”¹



With Provider 790 monitoring your fall risk patients, staff are free for more direct interaction with all of their patients.

The Provider 790 nurse call system's fall prevention features let staff remain completely mobile while continuously monitoring 'fall risk' patients and receiving alerts to critical events such as an at risk patient getting out of bed. As part of a patient safety program, **Provider 790 significantly increases patient safety through a combination of the following monitoring and streamlined alerting options:**



Patient and Bed Monitoring


Provider 790 lets staff proactively monitor fall risk patients in several ways: the dome light outside a patient's room indicates a fall risk, the patient's calls are noted as "fall risk", and their bed is continuously monitored for bed height, bed exit on/off, bed height and bed rails up/down.² If the patient gets out of bed, all caregivers on the patient's team are instantly notified of the "bed exit" alarm.

Bathroom Features for Fall Prevention

A significant fall risk is going to and from the bathroom. To make this trip much safer, Provider 790 lets patients specifically request assistance to the bathroom with the "toilet" call button. When they are finished, they can press the "CALL" button on the bathroom station and communicate their needs directly to their caregiver. If there is an emergency, a pull on the station's cord will send emergency alerts to everyone on the patient's team and (optionally) to every caregiver in the unit.



Automated patient rounding means one less thing for caregivers to remember.



Reducing potential falls with:


- Bed status monitoring and bed exit alarming
- Group notification of bed exit
- Bathroom assist notification for staff assistance to/from bathroom
- Automated rounding for consistent staff visits

Automated Patient Rounding

When patients know that their caregivers will consistently check on them, they are much less likely to attempt tasks on their own which is when falls can occur. The Provider 790 rounding feature lets staff set up automatic rounding reminders customized on a per patient basis. Caregivers are automatically alerted when a rounding period has expired; if the first caregiver is busy then other team members are automatically notified.

Wireless Notification

The Provider 790's rounding, audio bathroom, and patient bed monitoring functionality is greatly enhanced when integrated with wireless smart phone, wireless SIP phone, or pocket pager alerting to each patient's assigned caregivers. Over their phone or pager, caregivers receive notifications from their patients such as a "bed exit" alarm or an expired rounding reminder from anywhere within the facility – speeding notification and in turn the caregiver's response. From the bedside and bathroom intercom stations, patients tell their caregiver exactly what they need over the caregiver's wireless phone. The Provider 790 also includes group notification options to route high priority calls, such as "bed exit" or "bath emergency" to multiple caregiver's wireless devices so the closest caregiver can quickly respond.



To see how Provider 790 can be an integral component of your patient safety program, call **Jeron** at **1-800-621-1903** or visit us at **www.jeron.com**