



Explore
the
Options



SIP Wireless
Phones



Cordless
Phones



Text Messaging
& Pocket Paging



Real Time
Locating



PC Console
Display



Activity Logging
& Reporting



Automated
Overhead
Paging

NURSE STATION OPTIONS



Touchscreen Nurse Console

- Graphic display of multiple calls
- Embedded technology designed for continuous operation
- Wall or desk-mount option
- Menu-driven touchscreen navigation
- One-touch operation to alert staff
- Communicate with associated stations
- Tone-Visual option

Satellite Nurse Console

- Full functionality in a space saving design
- Wall or desk-mount option
- Multiple call-in display



Annunciator Panel

- Dedicated call indications for each room
- Full English call display
- Wall or desk-mount option



STAFF AREA OPTIONS

Duty Stations

- Tone/visual indication of outstanding calls
- Tone mute (regenerates with next call)
- Intercom option (to Nurse Consoles)



PATIENT/RESIDENT ROOM STATIONS

Patient Stations

- Built-in LED night light
- Normal, Priority, Personal Attention, and Cord Out call priorities
- Ultra-reliable tilt release connector
- Compatible with most bed siderail controls
- Entertainment and light control interface
- Tone-visual station options
- No dummy plug required



Pillow Speakers and Call Cords

- Completely sealed – easily cleaned
- UL listed for use in oxygen-enriched atmosphere
- Large easy-to-press buttons



Bath and Special Function Stations

- Built-in LED night light
- One-touch/pull operation
- Auxiliary alarm input (vent alarm, bed exit, IV alarm)
- Audio to the bathroom option

The Jeron Provider 680 Nurse Call System with system integrations **is a complete patient/resident to staff alerting and communications solution** for a nursing unit, a single building facility, or a multi-building campus. The integrated Provider 680 solution ensures direct alerting and communication for a safer environment, that drives both patient and resident satisfaction.

Jeron and our family of local factory certified and trained distributors provide a complete turn-key integrated solution for your facility with timely local service and support.

To find out more about the Jeron Provider 680 Nurse Call System, call **Jeron at 1-800-621-1903** or visit us at **www.jeron.com**