

Streamlined Notification to Speed Response

Multiple Modes of Annunciation Alert Caregivers to Critical Events

When a patient is in duress, code blue and rapid response alarms will mobilize a team of caregivers and direct them immediately to the patient. Provider 700's automated facility-wide alerting options can speed response times by replacing the traditional manual methods of calling the phone operator to dispatch crash teams.

Centralized Notification

Supporting instant code call and rapid response notification at the facility's switchboard or centralized dispatch area, Provider 700 displays these urgent events at one or more locations through networked large screens or centralized staff console. The code and rapid response tones immediately direct attention to the active event while large screens display the room number and descriptive location. From central staff consoles, urgent calls from anywhere throughout the facility can be answered and confirmed before dispatching staff.

When Seconds Count... Provider 700 Gives You the Tools to Respond Quickly!

Speed Response With:

- Centralized notification
- Automated overhead pages
- Wireless crash team alerts



Facility-Wide Code &
Rapid Response Display
at a Central Location

Wireless Crash Team Notification

Receiving critical information can be delayed or missed if your facility relies solely on public address announcements for emergency events. Provider 700's wireless alerts go directly to caregivers' smartphones, in-house SIP wireless phones, and/or pocket pagers. This means everyone on the crash team gets the same information at the same time. A single button press from anywhere in the facility instantly alerts the crash team with the room number and bed over their wireless devices. Provider 700's unique Group Alert technology alerts multiple smart or SIP phones simultaneously to an emergency event, eliminating the need for expensive middleware for group alerts.

These alerting options can be extended to other types of urgent patient and staff calls such as unit-wide "Staff Emergency" and "Bed Exit" alerts. Creating any number of custom team alerts ensures caregivers are immediately notified of all critical events so they can respond as quickly as possible.



Crash Team Stations

- Everywhere-access, large Code button stations within easy reach for an emergency
- Customizable call buttons to meet any type of team notification
- Distinct Dome Light indicators to direct the crash team



Automated Announcements

Precious seconds can be saved when crash team public address announcement are automated using Provider 700's automated paging feature. Provider 700 integrates with your facility's public address system for instant announcements of crash team events including the room number, nursing unit name, location within facility, and instructions. Multiple zone outputs give you the flexibility to determine what areas should receive the announcement, while the "call canceled" announcement lets mobile caregivers know that an accidental call has been cleared.



Learn more about the Provider 700 Nurse Call System; contact Jeron at 800-621-1903 or www.jeron.com